



2483 47th St.
Astoria, NY 11103
contact@bartovation.com

Buffer Tablets pH 4, 7, & 10 COC

LPH4710K30

For water analysis and calibration of pH meters in laboratories

These tablets come in sleeves of ten tablets of the same pH value and have a simple preparation process.

These tablets are manufactured under tightly controlled conditions and on state of the art machinery. These tightly controlled conditions in addition the hermetically sealed packaging of the tablets enables these tablets to have a high degree of precision and accuracy – for three to seven years from the date of manufacturing. For use, simply push the tablet from the foil directly into the sample, eliminating the need for the user to touch the buffer and ensuring the entire tablet enters the calibration solution.

For usage, simply dissolve one tablet into 20 mL of DI water. When done according to the instructions, the sample will be:

4: 4.0 +/- 0.05 pH; pink

7: 7.0 +/- 0.05 pH; yellow

10: 10.0 +/- 0.05 pH; green

All products released for sale are in conformance with the stated parameters and tested for precision and accuracy.

Scott Bartman 03/04/2020

Scott Bartman | Owner

Bartovation LLC.