



Bartovation LLC.
2483 47th St.
Astoria NY 11103
contact@bartovation.com

Product Specifications

Product: A test kit for the determination of somatic cell counts in fresh milk.

Test Strip Dimensions: Approximately. 25" width x 3" length

Activator Solution: Buffered solution containing preservative in a plastic dropper bottle and plastic cap.

Pipettes: 50 x Exact Volume Pipettes

Color Chart: Printed color chart for estimating results using color comparison.

Package Insert: instructions

Shelf Life: 12 months from date of manufacture

Sample Volume: 60 μ L milk

Activator Solution: After the milk is fully absorbed, apply 2 drops of buffer solution to the test pad.

Reaction Time: 15 minutes

Analytical Range: Color Chart values at 100, 200, 400, 800 (x1000) cells/ml